### Achtung

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### Caution

This product is restricted to use by or on the order of expert personnel / physicians.

### Symbols

 CE-sign and identification number of the Notified Body.  
The product is in conformity with the basic requirements of the Regulation 2017/745.  
**SN** Serial number

### Intended Use

The WISAP heating tube 7642HSE (incl. 7642HSE1, 76HSE10) in connection with the system insufflator and the Flow Therme of the WISAP Medical Technology serves for building up a pneumoperitoneum for diagnostic or operative Pelviscopy/Laparoscopy. The insufflation gas is preheated to body temperature by the heating tube. The heating tube has to be electrically operated with original WISAP appliances only. See also Instructions for use of the devices mentioned above.

### Contraindications

Use of the heating tube is contraindicated, whenever Pelviscopy/Laparoscopy is contraindicated. The heating tube 7642HSE (incl. 7642HSE1, 76HSE10) is only intended for the use with the above named devices

### Important Hints

- Use heating tube 7642HSE (incl. 7642HSE1, 76HSE10) for its intended purpose only.
- The sterile disposable heating tube is supplied sterile and has been sterilized with ethylene oxide
- Before use, the packaging must be checked for damage. If there is visible damage of any kind, the heating tube must not be used!
- The disposable heating tube must be used immediately after the sterile packaging. After accidental opening, the heating tube must be disposed of!
- The heating tube has to be used for heating gas (CO<sub>2</sub>) only.  
**The heating tube is not suitable for warming up liquids.**
- Prior to the operative Intervention, check all accessories as to their proper function. If faults are found or assumed, their use is strictly prohibited!
- Always use a new sterile heating tube and a new filter for each patient.
- Interventions requiring the use of the products described in these instructions, may only be performed by persons with adequate training and experience.
- Prior to each Intervention, please consult the pertinent medical literature for techniques, complications and hazards.
- Independent from the general recommendations in this manual, the selection of accessories needed for the operative Intervention is within the sole responsibility of the operating surgeon.
- For avoiding cross-contamination, back-flow and contamination, the original WISAP filter (REF 7070TFW2) has to be used. Filter and heating tube have to be replaced for each patient.  
**The filter is disposable and thus not re-sterilisable!**
- Do not bend, squeeze or over-stretch heating tube and do not roll it up with a diameter smaller than 25 cm.

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**WE RESERVED THE RIGHT FOR CHANGES ON  
TECHNICAL FEATURES AND DOCUMENTATION!**

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## Cleaning, Disinfection and Sterilization

### ATTENTION!

The heating tube 7642 HSE (incl. 7642HSE1, 76HSE10) is a disposable product and must not be treated.


After using this disposal must be supplied.

## Description of heating tube connection

see also the drawing shown below and the Instructions for use of the devices System Insufflator and Flow Therme

1. Connection socket for the gas
2. Connection plug of the heating tube cable
3. Gas-output fitting to the patient at the insufflator
4. Filter (disposable)
5. Filter adaptor
6. Gas-Coupling fitting at the heating tube
7. Gas heating tube
8. Heating tube coupling to instrument
9. Veres needle

Electrical Connection of the Heating Tube: 19,5 V/1 A DC

## Disposal of the heating tube

After usage the heating tube must be disposed through the clinical disposal system.

## Warranty

See out General Terms and Conditions. In case of unauthorized opening / modifications and/or repairs of the heating tube we cannot be held liable for proper and safe function of the device. All warranty claims are declared null and void then.

## INCIDENTS

**All serious incidents related to the device must be reported to the manufacturer and the component authority of the Member State in which the user and/or the patient is established.**