

# C<sup>3</sup> ECO4

MOBILE THERMAL ABLATION DEVICE  
FOR LOW RESSOURCE SETTINGS

# WISAP<sup>®</sup>



TIMER FUNCTION  
with visual and acoustic  
feedback

TOGGLE HEATING  
for power saving mode

FULL SET WITH 4 PROBES,  
BATTERY AND CASE

C3 ECO4 SET  
Art.-No. 6400

## C<sup>3</sup> ECO4 ADVANTAGES

- ✓ No gas required
- ✓ Fully mobile
- ✓ No autoclave required
- ✓ Simple handling
- ✓ Timer function
- ✓ Heat protection
- ✓ Different therapy probes  
(flat/nipple)

## CLINICAL BACKGROUND

### Treatment of cervical intraepithelial neoplasia (CIN)

The most important cause for the development of CIN is a chronic infection with human papillomavirus (HPV). HPV is very common in the population, particularly among younger women. Regardless of screening strategy, access to effective treatment for women with cervical pre-cancer is crucial. In some countries, up to 80% of women diagnosed with cervical pre-cancer do not receive the recommended treatment.

Therefore it is crucial to implement a clinically tested and highly efficacious treatment method, which is easy to handle and cheap to run.

# C<sup>3</sup> ECO4 MOBILE SOLUTION

SAFE, FAST AND EFFICACIOUS PREVENTION  
OF CERVICAL CANCER



13485:  
2016



## THE MOBILE SOLUTION

The C3 thermo coagulator is optimized for use in low-resource settings. It is independent of power grid and clinical infrastructure, it is handheld, mobile and easy to use. This LMIC-adapted thermo coagulator is designed with one idea in mind: Simple handling for maximum effect. Therefore the device features a single optimally designed thermo-probe, together with a fixed temperature setting and a timer function. That way we eliminate any possible source for human error and ensure the greatest efficacy on site.

### TECHNICAL DATA

Mains connection:	100 - 240V AC. 50/60Hz
Power input:	70 VA
Treatment Temperature:	100°C +/-10%
Dimensions:	w x h x d = 30 x 280 x 300 mm
Material Handle:	Makrolon, white PC
Material Thermo-probe:	PVC-U
Material anti-stick-coating:	Teflon (FDA)
Weight handle:	0,2 kg
Weight thermo-probe:	0,1 kg



### ORDER DATA CERVICAL COLD COAGULATOR / C3 ECO4

C3 ECO4 set				
6400	●	●●	●	●