



## THE MOBILE SOLUTION

Safe, fast and efficacious prevention of cervical cancer

# WISAP®



C3 - Cervical Cold Coagulator  
Art.-No. 610x

## C3 ADVANTAGES

- ✓ No gas required
- ✓ No power-grid required
- ✓ No autoclave required
- ✓ Simple handling
- ✓ Timer function
- ✓ Heat protection

## CLINICAL BACKGROUND

### Treatment of cervical intraepithelial neoplasia (CIN)

The most important cause for the development of CIN is a chronic infection with human papillomavirus (HPV). HPV is very common in the population, particularly among younger women.


Regardless of screening strategy, access to effective treatment for women with cervical pre-cancer is crucial. In some countries, up to 80% of women diagnosed with cervical pre-cancer do not receive the recommended treatment.

Therefore it is crucial to implement a clinically tested and highly efficacious treatment method, which is easy to handle and cheap to run.



# THE MOBILE SOLUTION

Safe, fast and efficacious prevention of cervical cancer



## THE MOBILE AND STURDY SOLUTION

The C3 thermo coagulator is optimized for use in low-resource settings. It is independent of power grid and clinical infrastructure, it is handheld, mobile and easy to use. This LMIC-adapted thermo coagulator is designed with one idea in mind: Simple handling for maximum effect. Therefore the device features a single optimally designed thermo-probe, together with a fixed temperature setting and a timer function. That way we eliminate any possible source for human error and ensure the greatest efficacy on site.

TECHNICAL DATA	
Mains connection:	100 - 240V AC. 50/60Hz
Power input:	70 VA
Treatment Temperature:	100°C +/-10%
Dimensions:	w x h x d = 30 x 280 x 300 mm
Material Handle:	Makrolon, white PC
Material Thermo-probe:	PPSU / VA
Material anti-stick-coating:	Teflon (FDA)
Weight handle:	0,2 kg
Weight thermo-probe:	0,1 kg



### ORDER DATA

### CERVICAL COLD COAGULATOR / C3

C3 set	C3 handle	Therapy probe 20 flat	Therapy probe 20 nipple	Therapy probe 17 flat	Slider heat-protection short 20	Slider heat-protection short 17
	6110	6120	6120N	6122	6112	6114
6100	●	●●			●●	
6102	●		●●	●●		●●
6104	●	●	●	●	●	●
6105	●	●	●		●●	
6106	●		●	●	●	●

All sets include a power supply (Article-Nr: 6130) with adapter plug

