

WISAP®

SafeCut

SMART MORCELLATION
safe and efficient tissue removal



Morcellation system
for safe and efficient tissue
removal in myomectomy
and hysterectomy.



proven technology | efficient healing
WISAP® Quality made in Germany

WISAP®



THE SMART HYBRID CONCEPT performance meets sustainability

Sterile single-use cutting tubes.

Always sharp, sterile single-use cutting tubes guarantee peak cutting performance in every procedure. A premium reusable motor, validated for up to 500 sterilization cycles, dramatically reduces per-case costs and environmental waste. This hybrid approach offers the ideal balance of cutting efficiency, reliability and economic value.

Adjustable speed. Ultimate control.

With a fully adjustable RPM range from 200 to 500, surgeons can tailor the device's performance to every clinical scenario. Whether working delicately around critical structures or powering through dense tissue, this morcellator adapts seamlessly to your technique.

Unmatched torque. Power through all cases.

Experience the highest torque output in its class, designed to deliver unwavering power even in the most demanding myomectomy or hysterectomy cases.

One Footpedal. Total control.

The intuitive single footswitch interface provides effortless control, allowing surgeons to keep their hands — and their focus — on the operative field. Smooth, ergonomic, and distraction-free.

Sustainable performance. No compromise.

By combining long-lasting technology with sterile single-use components, this system significantly cuts down medical waste compared to fully disposable alternatives — while maintaining the highest surgical standards.

ORDER DATA SAFECUT CONTROL-UNIT (SET)

7690S	SafeCut control-unit (set) Hybrid morecellation system for in-bag morcellation	7690MS	SafeCut procedure set Single use, 5 pcs / pack
7690BU	SafeCut control-unit with touch-display	7690CT	SafeCut cutting-unit, handle, cutting tube 15, double sealing Single use, double packed
7690M	SafeCut motor, autoclavable	7690B	SafeCut morcellation-bag, 2 litre volume, 2 entry-ports Single use, single packed
7690C	SafeCut cabel, autoclavable	1110ND2	Mains cable, L = 2,5 m / Typ D
7690FS	SafeCut footswitch, electrical	76900B	SafeCut obturator, autoclavable

WISAP® Medical Technology GmbH
Fichtenstraße 27
D-85649 Brunthal/Hofolding

Phone +49(0)8104 8908-0
Fax +49(0)8104 8908-90
wisap.de • info@wisap.de