

WISAP®

MORCELLATION

POWER DRIVE (SET) / ACCESSORIES



POWER DRIVE (SET)



MULTI-USE CUTTING TUBES



WISAP POWER DRIVE (Set)

| PRODUCT OVERVIEW | FEATURES |
|-------------------------|----------|
| Cutting tubes (Ø in mm) | 15 / 20 |
| Direct tissue removal | ● |
| Multi-use cutting tubes | ● |


POWER DRIVE (SET) – CUTTING TUBES



THE WISAP® MORCELLATION-SET

"Morcellation" describes a special cutting technique whereby large amounts of tissue can be removed via a minimally invasive approach in a short time. This method is used for both diagnostic purposes as well as indications for surgical intervention.

The WISAP® Morcellation-Set includes the morcellator-unit, matching cutting-tubes, instruments and trocars. All system components are designed to fit each other. Thus we provide a professional, impeccable system.



POWER DRIVE (SET)

POWER DRIVE (Set)

- Gearbox fits all sizes of morcellator-blades (Ø15 / 20)
- Multi-use cutting tubes



MULTI-USE
CUTTING TUBES

Ø MORCELLATION TUBES
15/20

TECHNICAL DATA POWER DRIVE (SET)

| | |
|-------------------|---|
| Controlling unit | |
| Mains connection: | 100 - 240 VAC |
| Power input: | 110 VA |
| Dimensions: | w x h x d = 350 x 90 x 290 mm |
| Weight: | 2,8 kg |
| Casing: | Stainless steel, aluminium / white coated |
| Drive Motor | |
| Dimensions: | l x w = 175 x 41 mm |
| Diameter: | 34 mm |
| Weight: | 0,62 kg |



ORDER DATA POWER DRIVE (SET)




| | | |
|----------|--------------------|-------------------------|
| | POWER DRIVE (Set) | |
| 7688 PD1 | Set consisting of: | |
| | 7688 PDU | Controlling unit |
| | 7688 PDC | Instrument cable, L 3 m |
| | 7688 PDM0 | Drive motor |
| | 7688 PDG1 | Gearbox Macro |



POWER DRIVE (SET) CUTTING TUBES / TROCARS








ORDER DATA

CUTTING TUBES





| MULTI USE | Cutting tube | | Seal (single) / 10 pcs. p. pack | |
|-----------|---|--------|---------------------------------|-------|
| Ø in mm | Medical grade steel / autoclaveable | | Silicon / Autoclaveable | |
| 15 |  | 7689M1 | 7620DK1 | white |
| |  | 7689W1 | 7620DK1 | white |
| 20 |  | 7689M2 | 7620DK1 | white |
| |  | 7689W2 | 7620DK1 | white |

 serrated edge  waved edge

cutting tubes
Ø 15 mm

| |  |  |  |  |  |
|-----------|--|--|--|--|--|
| Ø 15 mm | Flap valve | Conical tip | Triangular tip | Sleeve polished | Sleeve with thread |
| | 7676KF30 | 7676T3N | 7676T4N | 7676R1N | 7676R3N |
| 7676K4NO | ● | | ● | ● | |
| 7676K3NO | ● | ● | | ● | |
| 7676K4NOG | ● | | ● | | ● |
| 7676K3NOG | ● | ● | | | ● |

cutting tubes
Ø 20 mm



| |  |  |  |  |
|-----------|---|---|---|---|
| Ø 20 mm | Flap valve | Dilatation set | Sleeve polished | Sleeve with thread |
| | 7676KF50 | 7670-7N | 7676R2N | 7676R4N |
| 7676K5DN | ● | ● | ● | |
| 7676K5DNG | ● | ● | | ● |

INSTRUMENTS FOR MORCELLATION



ORDER DATA

MYOMA DRILL




| | | |
|--|--|---------------------------|
|  mm |  mm | Myoma Drill, acc. to Semm |
| 420 mm | 5 mm | 7689 MB5 |





| | | |
|--|--|---------------------------|
|  mm |  mm | Myoma Drill, acc. to Semm |
| 420 mm | 10 mm | 7689 MB1 |

ORDER DATA


TENACULUM FORCEPS






| | | |
|--|--|---|
|  mm |  mm | Tenaculum Forceps with combination handle |
| 420 mm | 10 mm | 7689 T1 |



ORDER DATA

CLAW FORCEPS



| | | |
|--|--|--------------------------------------|
|  mm |  mm | Claw Forceps with combination handle |
| 420 mm | 10 mm | 7689 ZG1 |





| | | |
|--|--|-----------------------------------|
|  mm |  mm | Claw Forceps with scissors handle |
| 420 mm | 5 mm | 7689 ZG5 |

| | |
|------------------------------|------------------------|
| Acc. for 7689 ZG5, -ZG1, -T1 | Sealing 10 pcs.p.pack. |
| | 7652 DK |



INSTRUMENTS FOR MORCELLATION



2 PARTS

| ORDER DATA | | TENACULUM FORCEPS | |
|------------|--|---|-------------------|
| |  |  | |
| | Handle without ratchet | Handle with ratchet | Instrument insert |
| | 7301HN1 | 7301HNR1 | |



|  mm |  mm | Claw Forceps two rows of teeth, double-action | | |
|--|--|---|-----------|----------|
| 335 mm | 10 mm | 7370 ZG1N | 7370 ZG2N | 7370 ZGE |

WISAP Medical Technology GmbH

Fichtenstraße 27

D - 85649 Brunnthal/Hofolding

Tel. +49(0)81 04/8908-10

Fax. +49(0)81 04/8908-90

www.wisap.de · info@wisap.de

W i S A P[®]